

Product Datasheet

Rat IL10Ra ELISA Kit (orb437142)

Description	Rat Interleukin 10 Receptor Alpha (IL10Ra) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	G3V830
Sensitivity	0.063ng/mL
Entrez	117539
Expiration Date	Please enquire.