

Product Datasheet

Mouse IL11Ra ELISA Kit (orb437149)

Description	Mouse Interleukin 11 Receptor Alpha (IL11Ra) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q64385
Sensitivity	0.061ng/mL
Entrez	16157
Expiration Date	Please enquire.