


Product Datasheet

Cytokeratin 5/18 antibody (orb43715)

Description	Mouse monoclonal antibody to Cytokeratin 5/18
Reactivity	Bovine, Canine, Hamster, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, IP, WB
Immunogen	Cytoskeleton preparation from HeLa cells.
Target	Cytokeratin 5/8
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive control: human prostate.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive control: human prostate.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of human p...