



Product Datasheet

Cytokeratin 8 antibody (orb43716)

Description	Mouse Monoclonal to Cytokeratin 8.
Reactivity	Bovine, Human, Porcine, Sheep
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, IP, WB
Immunogen	Cytoskeleton preparation from HeLa cells.
Target	Cytokeratin 8
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q7L4M3
Entrez	3856
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blotting
detection of
cytokerati...



Immunocytochemistry
staining of cytoker...