

Product Datasheet

Human IL12Rb2 ELISA Kit (orb437162)

Description	Human Interleukin 12 Receptor Beta 2 (IL12Rb2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q99665
Sensitivity	0.055ng/mL
Entrez	3595
Expiration Date	Please enquire.