

Product Datasheet

Canine IL13 ELISA Kit (orb437164)

Description	Canine Interleukin 13 (IL13) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9N0W9
Sensitivity	5.5pg/mL
Entrez	442990
Expiration Date	Please enquire.