

Product Datasheet

Mouse IL13Ra2 ELISA Kit (orb437172)

Description	Mouse Interleukin 13 Receptor Alpha 2 (IL13Ra2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O88786
Sensitivity	0.126ng/mL
Entrez	16165
Expiration Date	Please enquire.