

Product Datasheet

Bovine IL15 ELISA Kit (orb437173)

Description	Bovine Interleukin 15 (IL15) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q28028
Sensitivity	6.0pg/mL
Entrez	281248
Expiration Date	Please enquire.