

Product Datasheet

Bovine IL16 ELISA Kit (orb437178)

Description	Bovine Interleukin 16 (IL16) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q0V8R5
Sensitivity	6.1pg/mL
Entrez	506314
Expiration Date	Please enquire.