

Product Datasheet

Chicken IL16 ELISA Kit (orb437179)

Description	Chicken Interleukin 16 (IL16) Sandwich ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q70XC4
Sensitivity	5.7pg/mL
Entrez	374270
Expiration Date	Please enquire.