

Product Datasheet

Mouse IL16 ELISA Kit (orb437181)

Description	Mouse Interleukin 16 (IL16) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O54824
Sensitivity	9.12pg/mL
Entrez	16170
Expiration Date	Please enquire.