

Product Datasheet

Rat IL17B ELISA Kit (orb437187)

Description	Rat Interleukin 17B (IL17B) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G3V8K9
Sensitivity	6.4pg/mL
Entrez	116472
Expiration Date	Please enquire.