

Product Datasheet

Canine IL17 ELISA Kit (orb437188)

Description	Canine Interleukin 17 (IL17) Sandwich ELISA Kit
Reactivity	Canine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	C6L8D7
Sensitivity	31pg/mL
Entrez	481837
Expiration Date	Please enquire.