

Product Datasheet

Mouse IL18BP ELISA Kit (orb437204)

Description	Mouse Interleukin 18 Binding Protein (IL18BP) Sandwich ELISA Kit
Reactivity	Mouse
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9Z0M9
Sensitivity	45pg/mL
Entrez	16068
Expiration Date	Please enquire.