

Product Datasheet

Chicken IL18 ELISA Kit (orb437207)

Description	Chicken Interleukin 18 (IL18) Sandwich ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8QFQ8
Sensitivity	7.3pg/mL
Entrez	395312
Expiration Date	Please enquire.