

Product Datasheet

Rabbit IL18 ELISA Kit (orb437216)

Description	Rabbit Interleukin 18 (IL18) Sandwich ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	B1A3U4
Sensitivity	6.2pg/mL
Entrez	100144338
Expiration Date	Please enquire.