

Product Datasheet

Equine IL1a ELISA Kit (orb437222)

Description	Equine Interleukin 1 Alpha (IL1a) Sandwich ELISA Kit
Reactivity	Equine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q28385
Sensitivity	27pg/mL
Entrez	100064969
Expiration Date	Please enquire.