

Product Datasheet

Guinea pig IL1a ELISA Kit (orb437224)

Description	Guinea pig Interleukin 1 Alpha (IL1a) Sandwich ELISA Kit
Reactivity	Guinea pig
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	H0W0J5
Sensitivity	6.6pg/mL
Entrez	100135555
Expiration Date	Please enquire.