

Product Datasheet

Rabbit IL1a ELISA Kit (orb437229)

Description	Rabbit Interleukin 1 Alpha (IL1a) Sandwich ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P04822
Sensitivity	5.4pg/mL
Entrez	100009250
Expiration Date	Please enquire.