

Product Datasheet

Bovine IL1 beta ELISA Kit (orb437230)

Description	Bovine IL1 beta ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P09428
Sensitivity	6.1pg/mL
Entrez	281251
Expiration Date	Please enquire.