

Product Datasheet

Chicken IL1 beta ELISA Kit (orb437232)

Description	Chicken IL1 beta ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	073909
Sensitivity	5.7pg/mL
Entrez	395196
Expiration Date	Please enquire.