

Product Datasheet

Guinea pig IL1 beta ELISA Kit (orb437234)

Description	Guinea pig IL1 beta ELISA Kit
Reactivity	Guinea pig
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9WVG1
Sensitivity	25pg/mL
Entrez	100135556
Expiration Date	Please enquire.