

Product Datasheet

Rat IL1 beta ELISA Kit (orb437238)

Description	Rat IL1 beta ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q63264
Sensitivity	5.5pg/mL
Entrez	24494
Expiration Date	Please enquire.