

Product Datasheet

Mouse IL1R2 ELISA Kit (orb437247)

Description	Mouse Interleukin 1 Receptor Type II (IL1R2) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P27931
Sensitivity	11.38pg/mL
Entrez	16178
Expiration Date	Please enquire.