

Product Datasheet

CD28 antibody (FITC) (orb43725)

Description

Mouse Monoclonal to CD28.

Reactivity

Human, Primate

Conjugation

FITC

Tested Applications

FC

Immunogen

DC28.1.3.3 murine T cell hybridoma transfected with human CD28 cDNA

Target

CD28

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P10747

Entrez


940

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.



Separation of human CD28 positive lympho...



Flow cytometry surface staining pattern ...