

Product Datasheet

Human IL1RAP ELISA Kit (orb437254)

Description	Human Interleukin 1 Receptor Accessory Protein (IL1RAP) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9NPH3
Sensitivity	0.065ng/mL
Entrez	3556
Expiration Date	Please enquire.