

Product Datasheet

Mouse IL1RAP ELISA Kit (orb437256)

Description	Mouse Interleukin 1 Receptor Accessory Protein (IL1RAP) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q61730
Sensitivity	14.2pg/mL
Entrez	16180
Expiration Date	Please enquire.