

Product Datasheet

Rat IL1RA ELISA Kit (orb437257)

Description	Rat Interleukin 1 Receptor Antagonist (IL1RA) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P25086
Sensitivity	12.1pg/mL
Entrez	60582
Expiration Date	Please enquire.