

Product Datasheet

Rat IL22Ra2 ELISA Kit (orb437273)

Description	Rat Interleukin 22 Receptor Alpha 2 (IL22Ra2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q7TNI4
Sensitivity	0.105ng/mL
Entrez	444986
Expiration Date	Please enquire.