

Product Datasheet

Human IL23R ELISA Kit (orb437277)

Description	Human Interleukin 23 Receptor (IL23R) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q5VWK5
Sensitivity	0.041ng/mL
Entrez	149233
Expiration Date	Please enquire.