

Product Datasheet

Rat IL23R ELISA Kit (orb437278)

Description	Rat Interleukin 23 Receptor (IL23R) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	F1LX96
Sensitivity	0.042ng/mL
Entrez	689932
Expiration Date	Please enquire.