

Product Datasheet

Canine IL2 ELISA Kit (orb437293)

Description	Canine Interleukin 2 (IL2) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q29416
Sensitivity	6.2pg/mL
Entrez	403989
Expiration Date	Please enquire.