

Product Datasheet

Chicken IL2 ELISA Kit (orb437294)

Description	Chicken Interleukin 2 (IL2) Sandwich ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O42288
Sensitivity	6.5pg/mL
Entrez	373958
Expiration Date	Please enquire.