

Product Datasheet

Equine IL2 ELISA Kit (orb437295)

Description	Equine Interleukin 2 (IL2) Sandwich ELISA Kit
Reactivity	Equine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P37997
Sensitivity	10.5pg/mL
Entrez	100034204
Expiration Date	Please enquire.