

Product Datasheet

Human IL2Rb ELISA Kit (orb437304)

Description	Human Interleukin 2 Receptor Beta (IL2Rb) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P14784
Sensitivity	28pg/mL
Entrez	3560
Expiration Date	Please enquire.