

Product Datasheet

Human IL2Rg ELISA Kit (orb437305)

Description	Human Interleukin 2 Receptor Gamma (IL2Rg) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P31785
Sensitivity	0.067ng/mL
Entrez	3561
Expiration Date	Please enquire.