

Product Datasheet

Mouse IL33 ELISA Kit (orb437311)

Description	Mouse Interleukin 33 (IL33) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8BVZ5
Sensitivity	31pg/mL
Entrez	77125
Expiration Date	Please enquire.