

Product Datasheet

Mouse IL34 ELISA Kit (orb437314)

Description	Mouse Interleukin 34 (IL34) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8R1R4
Sensitivity	5.3pg/mL
Entrez	76527
Expiration Date	Please enquire.