

Product Datasheet

Human IL4I1 ELISA Kit (orb437330)

Description	Human Interleukin 4 Induced Protein 1 (IL4I1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96RQ9
Sensitivity	0.052ng/mL
Entrez	259307
Expiration Date	Please enquire.