

Product Datasheet

Human IL4R ELISA Kit (orb437335)

Description	Human Interleukin 4 Receptor (IL4R) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P24394
Sensitivity	5.8pg/mL
Entrez	3566
Expiration Date	Please enquire.