

Product Datasheet

Porcine IL5 ELISA Kit (orb437345)

Description	Porcine Interleukin 5 (IL5) Sandwich ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9MYM5
Sensitivity	6.5pg/mL
Entrez	397409
Expiration Date	Please enquire.