

Product Datasheet

Human IL5Ra ELISA Kit (orb437347)

Description	Human Interleukin 5 Receptor Alpha (IL5Ra) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q01344
Sensitivity	24.6pg/mL
Entrez	3568
Expiration Date	Please enquire.