

Product Datasheet

C-peptide of Proinsulin antibody (Biotin) (orb43735)

Description	Mouse Monoclonal to C-peptide of Proinsulin.
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-P, RIA
Immunogen	Human C-peptide conjugated to bovine serum albumin.
Target	C-peptide
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Expiration Date	12 months from date of receipt.