

Product Datasheet

Equine IL6 ELISA Kit (orb437352)

Description	Equine Interleukin 6 (IL6) Sandwich ELISA Kit
Reactivity	Equine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	F6X6K9
Sensitivity	28pg/mL
Entrez	100034196
Expiration Date	Please enquire.