

Product Datasheet

Goat IL6 ELISA Kit (orb437353)

Description	Goat Interleukin 6 (IL6) Sandwich ELISA Kit
Reactivity	Goat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	E3UT57
Sensitivity	6.1pg/mL
Entrez	100860785
Expiration Date	Please enquire.