

Product Datasheet

Human IL6R ELISA Kit (orb437359)

Description	Human Interleukin 6 Receptor (IL6R) Sandwich ELISA Kit
Reactivity	Human
Range	1.875-120ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, urine and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, urine and other biological fluids.
Uniprot ID	P08887
Sensitivity	0.67ng/mL
Entrez	3570
Expiration Date	Please enquire.