

Product Datasheet

Sheep IL6 ELISA Kit (orb437362)

Description	Sheep Interleukin 6 (IL6) Sandwich ELISA Kit
Reactivity	Sheep
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P29455
Sensitivity	6.8pg/mL
Entrez	443406
Expiration Date	Please enquire.