

Product Datasheet

Rabbit IL8Ra ELISA Kit (orb437381)

Description	Rabbit Interleukin 8 Receptor Alpha (IL8Ra) Sandwich ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P21109
Sensitivity	0.061ng/mL
Entrez	100328622
Expiration Date	Please enquire.