

Product Datasheet

Mouse IL8Rb ELISA Kit (orb437384)

Description	Mouse Interleukin 8 Receptor Beta (IL8Rb) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35343
Sensitivity	12.4pg/mL
Entrez	12765
Expiration Date	Please enquire.