

Product Datasheet

Rat IL8Rb ELISA Kit (orb437385)

Description	Rat Interleukin 8 Receptor Beta (IL8Rb) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35407
Sensitivity	0.109ng/mL
Entrez	29385
Expiration Date	Please enquire.